

## Product Datasheet

# Human 6-O-Methylguanine-DNA Methyltransferase protein (orb669570)

<b>Description</b>	Recombinant Human Methylated-DNA--protein-cysteine methyltransferase is produced by our E.coli expre...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 20mMTris-HCl, 1mM EDTA, pH 8.0 .
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	N-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">P16455</a>
<b>Expiration Date</b>	6 months from date of receipt.