

Product Datasheet

Human Pro-Neuregulin-1 protein (orb669579)

Description	Recombinant Human Neuregulin-1 beta is produced by our E.coli expression system and the target gene ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 4mM HCl.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Uniprot ID	AAA58639.1
Expiration Date	6 months from date of receipt.