

## Product Datasheet

### Human CLDN18.2 protein (orb669651)

<b>Description</b>	Recombinant Human Claudin-18 Isoform A2 is produced by our E.coli expression system and the target g...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 25mM Tris-HCl, 25mM NaCl, 0.1% Triton X-100, 10% glycerol, pH 8.0 .
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	N-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P56856</b>
<b>Expiration Date</b>	6 months from date of receipt.