

Product Datasheet

DCIR antibody (orb669738)

Description	Mouse monoclonal antibody to DCIR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	This antibody was raised by immunising BALB/c mice with a dendritic cell immunoreceptor ectodomain-human Fc fusion protein.
Target	DCIR
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody has been used in flow cytometric analysis of DCIR expression by peripheral blood mononuclear cells (Klechevsky et al, 2010). A recombinant humanised form of this antibody fused to protein antigens is under investigation as an antigen-targeting vaccine (Flamar et al, 2012), and has been shown to trigger antigen-specific CD8+ T cell immunity.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q9UMR7
Expiration Date	12 months from date of receipt.