


Product Datasheet

LDLR antibody (orb66977)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal to LDLR.The low density lipoprotein receptor (LDLR) gene family consists of cell s... |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | N-terminal region of human LDLR |
| Preservatives | 0.01% NaN3 |
| Storage | -20°C |
| Note | For research use only |
| Clonality | Polyclonal |
| Purity | Affinity purification |
| MW | 95 |
| Uniprot ID | P01130 |
| Dilution Range | WB:1:1000-1:5000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of total cell extr...