

## Product Datasheet

### Ataxin 1 (Phospho-S775) antibody (orb669851)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to Ataxin 1 (Phospho-S775)   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IF, IHC, WB   |
| <b>Immunogen</b>           | Synthetic phosphopeptide derived from human Ataxin 1 around the phosphorylation site of Serine 775.                                     |
| <b>Preservatives</b>       | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2   |
| <b>Form/Appearance</b>     | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2   |
| <b>Concentration</b>       | 1 mg/ml   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | ~ 105 kDa   |
| <b>Uniprot ID</b>          | <a href="#">P54253</a>  |
| <b>Dilution Range</b>      | WB: 1:500-1:1000, IHC-P: 1:50-1:200   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |