

Product Datasheet

Histone H4 (Phospho-S47) antibody (orb669865)

Description	Rabbit polyclonal antibody to Histone H4 (Phospho-S47)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic phosphopeptide derived from human Histone H4 around the phosphorylation site of Serine 47.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 11 kDa
Uniprot ID	P62805
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.