

Product Datasheet

CDK5 (Phospho-Y15) antibody (orb669874)

Description	Rabbit polyclonal antibody to CDK5 (Phospho-Y15)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic phosphopeptide derived from human CDK5 around the phosphorylation site of Tyrosine 15.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 33 kDa
Uniprot ID	Q00535
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.