

Product Datasheet

MKK3 (Phospho-T222) antibody (orb669890)

Description	Rabbit polyclonal antibody to MKK3 (Phospho-T222)
Species/Host	Rabbit
Reactivity	Human, Mouse, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic phosphopeptide derived from human MKK3 around the phosphorylation site of Threonine 222.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 40 kDa
Uniprot ID	P46734
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.