

Product Datasheet

p47 phox (Phospho-S345) antibody (orb669908)

Description	Rabbit polyclonal antibody to p47 phox (Phospho-S345)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic phosphopeptide derived from human p47 phox around the phosphorylation site of Serine 345.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 47 kDa
Uniprot ID	P14598
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.