

Product Datasheet

GPR196 antibody (orb670178)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to GPR196 |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Synthetic peptide, corresponding to Human RGPL1/2/3/4/5/8. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.4. |
| Concentration | 1 mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 198 kDa |
| Uniprot ID | P0DJD0 , P0DJD1 , O14715 , Q99666 , A6NKT7 , Q7Z3J3 |
| Dilution Range | WB: 1:500-1000 |
| Expiration Date | 12 months from date of receipt. |