

Product Datasheet

Goat Anti-Rabbit IgG (H+L) Unconjugated (orb670230)

Description	Affinity purified goat anti-rabbit IgG(H & L)
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Rabbit IgG (H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, containing 1% rAlbumin
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Dilution Range	EIA and Western blots: 1:5,000-1:100,000. Immunohistochemistry: 1:500-1:5,000. Optimal working dilutions must be determined by end user.
Expiration Date	12 months from date of receipt.