

Product Datasheet

Goat Anti-Rabbit IgG (H+L) HRP (orb670232)

Description	Peroxidase-conjugated AffiniPure Goat Anti-rabbit IgG (H+L)
Species/Host	Goat
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Rabbit IgG (H+L)
Preservatives	0.02% thimerosal. 1mg/ml in PBS with 0.02% thimerosal, 50% glycerol, pH 7.4
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Dilution Range	EIA and Western blots: 1:5,000-1:100,000. Optimal working dilutions must be determined by end user.
Expiration Date	12 months from date of receipt.