

Product Datasheet

Goat Anti-Rabbit IgG (H+L) FITC (orb670237)

Description	Fluorescein-5-isothiocyanate(FITC-isomer 1) Amax=492nm;Emax=520nm Fluorophore/Protein: 8.5ug/mg;3.2 ...
Species/Host	Goat
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Rabbit IgG (H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6,containing 1% rAlbumin,50% glycerol
Concentration	1 mg/ml
Storage	Store at 4°C. Avoid freeze-thaw cycles. Prepare working dilution fresh each day.
Note	For research use only
Dilution Range	1:50 -1:200 for most applications
Expiration Date	12 months from date of receipt.