

Product Datasheet

Mouse Anti-Goat IgG (H+L) FITC (orb670238)

| | |
|----------------------------|--|
| Description | Fluorophore: Fluorescein Isothiocyanate (FITC) Fluorophore/Protein : 7.1 µg/mg; 2.7 moles FITC per ... |
| Species/Host | Mouse |
| Conjugation | FITC |
| Tested Applications | FC, IF |
| Immunogen | Goat IgG(H+L) |
| Preservatives | 0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, containing 1% rAlbumin, 50% glycerol |
| Concentration | 1.4 mg/ml |
| Storage | Store at 4°C. Avoid freeze-thaw cycles. Prepare working dilution fresh each day. |
| Note | For research use only |
| Dilution Range | 1:50 -1:200 for most applications |
| Expiration Date | 12 months from date of receipt. |