

## Product Datasheet

# Goat Anti-Rabbit IgG (H+L) Rhodamine (TRITC) (orb670240)

<b>Description</b>	Fluorophore: Tetramethyl rhodamine 5(and6)-isothiocyanate(TRITC) Amax=552nm;Emax=570nm Fluorophore/P...
<b>Species/Host</b>	Goat
<b>Conjugation</b>	TRITC
<b>Tested Applications</b>	FC, IF
<b>Immunogen</b>	Rabbit IgG (H+L)
<b>Preservatives</b>	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6,containing 1% rAlbumin,50% glycerol
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C. Avoid freeze-thaw cycles. Prepare working dilution fresh each day.
<b>Note</b>	For research use only
<b>Dilution Range</b>	1:50 -1:200 for most applications
<b>Expiration Date</b>	12 months from date of receipt.