

Product Datasheet

Goat Anti-Rabbit IgG (H+L) Biotin (orb670245)

Description	Biotin-SP-conjugated AffiniPure Goat Anti-Rabbit IgG (H+L)
Species/Host	Goat
Conjugation	Biotin
Tested Applications	ELISA, FC, ICC, IHC, WB
Immunogen	Rabbit IgG (H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, containing 1% rAlbumin, 50% glycerol
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Dilution Range	1:20000-1:400,000 for ELISA and Western blotting using enzyme-conjugated streptavidin; 1:500-1:5,000 for enzyme immunohisto/cytochemistry; 1:200-1:1,000 for flow cytometry and fluorescence immunohisto/cytochemistry
Expiration Date	12 months from date of receipt.