

Product Datasheet

Mouse Anti-Goat IgG (H+L) AP (orb670247)

Description	Alkaline Phosphatase-conjugated AffiniPure Mouse Anti-Goat IgG (H+L)
Species/Host	Mouse
Conjugation	AP
Tested Applications	ELISA, WB
Immunogen	Goat IgG(H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, containing 1% rAlbumin, 50% glycerol
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Dilution Range	EIA and Western blots: 1:5,000-1:100,000. Optimal working dilutions must be determined by end user.
Expiration Date	12 months from date of receipt.