

Product Datasheet

Goat Anti-Rabbit IgG (H+L) AP (orb670248)

Description	Alkaline Phosphatase-conjugated AffiniPure Goat Anti-Rabbit IgG (H+L)
Species/Host	Goat
Conjugation	AP
Tested Applications	ELISA, WB
Immunogen	Rabbit IgG (H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, containing 1% rAlbumin, 50% glycerol
Concentration	1 mg/ml
Storage	Store at 4°C. Avoid freeze-thaw cycles. Prepare working dilution fresh each day.
Note	For research use only
Dilution Range	1:5,000-1:50,000 for ELISA and Western blotting
Expiration Date	12 months from date of receipt.