

Product Datasheet

Rabbit Anti-Sheep IgG (H+L) FITC (orb670251)

Description	Fluorophore: Fluorescein Isothiocyanate (FITC) Fluorophore/Protein : 8.4 µg/mg; 3.2moles FITC per m...
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Sheep IgG (H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6,containing 1% rAlbumin,50% glycerol
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Dilution Range	1:50 -1:200 for most applications
Expiration Date	12 months from date of receipt.