

## Product Datasheet

### Goat Anti-Rat IgG(H+L) FITC (orb670252)

<b>Description</b>	Fluorophore: Fluorescein-5- isothiocyanate(FITC-isomer 1) Amax=492nm;Emax=520nm Fluorophore/Protein: ...
<b>Species/Host</b>	Goat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC, IF
<b>Immunogen</b>	Rat IgG (H+L)
<b>Preservatives</b>	0.05% Sodium Azide. 0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH 7.6
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Dilution Range</b>	1:50 -1:200 for most applications
<b>Expiration Date</b>	12 months from date of receipt.