

Product Datasheet

Goat Anti-Mouse IgG(H+L) Cy3 (orb670255)

Description	Fluorophore: Cy3 bis-NHS ester Amax=550 nm; Emax=570 nm Fluorophore/Protein: A552/A280=2.4
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Mouse IgG (H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH 7.6
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Dilution Range	1:100 - 1:800 for most applications
Expiration Date	12 months from date of receipt.