

Product Datasheet

Goat Anti-Rabbit IgG(H+L) Cy3 (orb670256)

| | |
|----------------------------|--|
| Description | Fluorophore: Cy3 bis-NHS ester Amax=550 nm; Emax=570 nm Fluorophore/Protein: A552/A280=2.2 |
| Species/Host | Goat |
| Conjugation | Cy3 |
| Tested Applications | FC, IF |
| Immunogen | Rabbit IgG (H+L) |
| Preservatives | 0.05% Sodium Azide. 0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH 7.6 |
| Concentration | 1 mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Dilution Range | 1:100 - 1:800 for most applications |
| Expiration Date | 12 months from date of receipt. |