

## Product Datasheet

### Goat Anti-Rat IgG(H+L) HRP (orb670289)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Peroxidase(HRP) -conjugated AffiniPure Goat Anti-Rat IgG (H+L)   |
| <b>Species/Host</b>        | Goat   |
| <b>Conjugation</b>         | HRP  |
| <b>Tested Applications</b> | ELISA, IHC, WB   |
| <b>Immunogen</b>           | Rat IgG (H+L)  |
| <b>Preservatives</b>       | None (Warning: Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.).<br>0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH 7.6 |
| <b>Concentration</b>       | 1 mg/ml  |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.   |
| <b>Note</b>                | For research use only  |
| <b>Dilution Range</b>      | 1:500 - 1:5000 for enzyme immunohisto/cytochemistry<br>1:5000 - 1:100000 for ELISA and Western blotting  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |