



Product Datasheet

RBM47 Antibody (orb670325)


Description	RBM47 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	E.coli-derived human RBM47 recombinant protein (Position: H66-K317).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64099 MW
Uniprot ID	A0AV96
Expiration Date	12 months from date of receipt.



Western blot analysis of RBM47 using ant...



IF analysis of RBM47 using anti-RBM47 an...



Flow Cytometry analysis of A431 cells us...