


Product Datasheet

BHLHE40 antibody (orb67045)

Description	Rabbit polyclonal to BHLHE40. BHLHE40 is a basic helix-loop-helix protein expressed in various tissue...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse BHLHE40
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	45
Uniprot ID	O35185
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...