


Product Datasheet

AGO3 antibody (orb67056)

Description	Rabbit polyclonal to AGO3. Protein argonaute-3 (AGO3) is also known as Eukaryotic translation initiat...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human AGO3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	97
Uniprot ID	Q9H9G7
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...