


Product Datasheet

EHMT2 antibody (orb67069)

Description	Rabbit polyclonal to EHMT2. Euchromatic histone-lysine N-methyltransferase 2 (EHMT2) belongs to the h...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human EHMT2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	132
Uniprot ID	Q96KQ7
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...