


Product Datasheet

MC4R antibody (orb67093)

Description	Rabbit polyclonal to MC4R.Melanocortin receptor 4 (MC4R) is a multi-pass membrane protein that belon...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MC4R
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	37
Uniprot ID	P32245
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...