


## Product Datasheet

### AXIN2 antibody (orb67132)

<b>Description</b>	Rabbit polyclonal to AXIN2. Axin-2 (AXIN2) may act as an inhibitor of the Wnt signaling pathway by do...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	C-terminal region of human AXIN2
<b>Preservatives</b>	0.01% NaN <sub>3</sub>
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>MW</b>	94
<b>Uniprot ID</b>	<a href="#">Q9Y2T1</a>
<b>Dilution Range</b>	WB:1:500-1:1000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of total cell extr...