


## Product Datasheet

### ATP11A antibody (orb67158)

<b>Description</b>	Rabbit polyclonal to ATP11A. Probable phospholipid-transporting ATPase 1H (ATP11A) is a multi-membran...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	middle region of human ATP11A
<b>Preservatives</b>	0.01% NaN <sub>3</sub>
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>MW</b>	130
<b>Uniprot ID</b>	<b>P98196</b>
<b>Dilution Range</b>	WB:1:500-1:2000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of total cell extr...