


Product Datasheet

MAN1a1 antibody (orb67191)

Description	Rabbit polyclonal to MAN1a1. Mannosyl-oligosaccharide 1,2- α -mannosidase IA (MAN1A1) is a type II ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MAN1a1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	73
Uniprot ID	P33908
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...