


Product Datasheet

MAN1a1 antibody pair (orb67193)

Description	Rabbit polyclonal to MAN1a1. 50 μl MAN1a1 (N-term1) &50 μl MAN1a1 (N-term2) .vPair antibodies ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	73
Uniprot ID	P33908
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...