


Product Datasheet

MMP17 antibody (orb67217)

Description	Rabbit polyclonal to MMP17.Matrix metalloproteinase-17 (MMP17) belongs to the peptidase M10A family....
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human MMP17
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	67
Uniprot ID	Q9ULZ9
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...