


Product Datasheet

MKRN1 antibody (orb67218)

Description	Rabbit polyclonal to MKRN1.E3 ubiquitin-protein ligase makorin-1 (MKRN1) catalyzes the covalent atta...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MKRN1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	53
Uniprot ID	Q9UHC7
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...