


Product Datasheet

MOCOS antibody (orb67220)

Description	Rabbit polyclonal to MOCOS. Molybdenum cofactor sulfurase (MOCOS) belongs to the class-V pyridoxal-ph...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MOCOS
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	98
Uniprot ID	Q96EN8
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...