



Product Datasheet

TET1 antibody pair (orb67233)

Description	Rabbit polyclonal to TET1. 50 μl TET1 (N-term1) &50 μl TET1 (C-term2) .vPair antibodies repres...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	235
Uniprot ID	Q8NFU7
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...



Western blot analysis of total cell extr...