


Product Datasheet

MST1R antibody (orb672377)

| | |
|----------------------------|---|
| Description | MST1R antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Synthesized peptide derived from the N-terminal region of Human Ron α . |
| Target | MST1R |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q04912 |
| Dilution Range | WB: 1:500-1:2000, ELISA: 1:10000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HeLa cells using...