


Product Datasheet

MLK1 antibody (orb67265)

Description	Rabbit polyclonal to MLK1. Mixed lineage kinase 1 (MLK1) is also known as Mitogen-activated protein k...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human MLK1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	122
Uniprot ID	P80192
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...