


Product Datasheet

IKKA antibody (orb67269)

Description	Rabbit polyclonal to IKKA. Inhibitor of nuclear factor kappa-B kinase subunit alpha (IKKA) belongs to...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human IKKA
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	85
Uniprot ID	O15111
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...