


Product Datasheet

PIK3AP1 antibody (orb67286)

Description	Rabbit polyclonal to PIK3AP1. Phosphoinositide 3-kinase adapter protein 1 (PIK3AP1) is involved in th...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal region of human PIK3AP1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	90
Uniprot ID	Q6ZUJ8
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...