


Product Datasheet

PGC antibody (orb673062)

Description	PGC antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PEP1
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04711
Dilution Range	WB: 1: 2000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Maize using PGC...