

## Product Datasheet

### CTU1 antibody (orb674946)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | CTU1 antibody  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, WB  |
| <b>Immunogen</b>           | Synthesized peptide derived from the C-terminal region of Human ATPBD3.  |
| <b>Target</b>              | CTU1   |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.   |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks.<br>For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">Q7Z7A3</a>   |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, ELISA: 1:20000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |