



Product Datasheet

CACNA1D antibody (orb676206)

Description	CACNA1D antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic Peptide
Target	CACNA1D
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01668
Dilution Range	IHC-P: 1: 100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse Br...



Immunohistochemical staining of Rat Brai...