


Product Datasheet

KCNA5 antibody (orb676229)

Description	KCNA5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide
Target	KCNA5
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22460
Dilution Range	WB: 1: 1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T using KCNA...