



## Product Datasheet

### IL31 antibody (orb676350)


|                            |                                                                              |
|----------------------------|------------------------------------------------------------------------------|
| <b>Description</b>         | IL31 antibody                                                                |
| <b>Species/Host</b>        | Rabbit                                                                       |
| <b>Reactivity</b>          | Human                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                 |
| <b>Tested Applications</b> | ELISA, WB                                                                    |
| <b>Immunogen</b>           | Synthetic peptide from Human protein at AA range: 21-70                      |
| <b>Target</b>              | IL31                                                                         |
| <b>Preservatives</b>       | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| <b>Form/Appearance</b>     | Liquid                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Note</b>                | For research use only                                                        |
| <b>Isotype</b>             | IgG                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q6EBC2</a>                                                       |
| <b>Dilution Range</b>      | WB: 1: 500-2000, ELISA: 1: 10000-20000                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                              |



Immunohistochemical staining of human-sp...



Western blot analysis of Hela Cell lysat...



Immunohistochemical staining of human-te...