



Product Datasheet

DNAJB1 antibody (orb676499)


Description	DNAJB1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Protein
Target	DNAJB1
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25685
Dilution Range	WB: 1:500-1:2000, IHC-P: 1: 100-1: 300, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Br...



Western blot analysis of 1) C2C12 Cell L...



Immunohistochemical staining of human Co...