

Product Datasheet

FNBP1L antibody (orb678362)

Description	FNBP1L antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Human FNBP1L
Target	FNBP1L
Preservatives	0.1M NaHCO ₃ , 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5TON5
Expiration Date	12 months from date of receipt.