

Product Datasheet

HECTD1 antibody (orb678735)

Description	HECTD1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Human HECTD1
Target	HECTD1
Preservatives	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C , Avoid freeze / thaw cycles.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9ULT8
Expiration Date	12 months from date of receipt.