

## Product Datasheet

### PRPH Antibody (orb680533)

|                      |                                                                                                                                                                                                       |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number       | orb680533                                                                                                                                                                                             |
| Category             | Antibodies                                                                                                                                                                                            |
| Description          | This PRPH Antibody is an unconjugated polyclonal product. It targets PRPH using a Human PRPH as the immunogen. This antibody is suitable for ELISA, IHC, WB. Purification: Antigen Affinity Purified. |
| Target               | PRPH                                                                                                                                                                                                  |
| Clonality            | Polyclonal                                                                                                                                                                                            |
| Species/Host         | Rabbit                                                                                                                                                                                                |
| Isotype              | IgG                                                                                                                                                                                                   |
| Conjugation          | Unconjugated                                                                                                                                                                                          |
| Reactivity           | Human, Mouse, Rat                                                                                                                                                                                     |
| Form/Appearance      | Liquid                                                                                                                                                                                                |
| Buffer/Preservatives | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.                                                                                                                  |
| Purification         | Antigen Affinity Purified                                                                                                                                                                             |
| Immunogen            | Human PRPH                                                                                                                                                                                            |
| UniProt ID           | <b>P41219</b>                                                                                                                                                                                         |
| Tested applications  | ELISA, IHC, WB                                                                                                                                                                                        |
| Storage              | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.                                                               |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Note**

For research use only

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)